



Product specification

Tube, 5 ml, (LxØ): 75 x 12 mm, PS



Product information

Product number	55.1579
----------------	---------

Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 12 mm, PS
Product number	55.1579
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, suitable for flow cytometer like FACSCanto, FACS Aria, FACSCalibur, LRS2, FC500, 500 piece(s)/StackPack

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	12 mm
Length of product	75 mm

Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	4000 x g

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	112000
Height of inner box	179 mm
Width of inner box	77 mm
Depth of inner box	394 mm
Height of case	175 mm
Width of case	383 mm
Depth of case	428 mm
Volume of case	0.0287 cbm
Weight of case	5.82 kg
EAN inner box	4038917157735
EAN case	4038917157728

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT