



Product specification

Tube, 13 ml, (LxØ): 95 x 16.8 mm, PS



Product information

Product number	55.468.001
----------------	------------

Product descriptions

Name	Tube, 13 ml, (LxØ): 95 x 16.8 mm, PS
Product number	55.468.001
Product description	Tube, working volume: 13 ml, (LxØ): 95 x 16.8 mm, material: PS, round base, transparent, push cap, natural, cap assembled, sterile, 100 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent
Color of cap	natural

Size

Diameter	16.8 mm
Length of product	95 mm

Product characteristics

Base shape	Round base
Graduation	no
Cap	Standard cap
Closure type	Push cap

Centrifugation max (RCF)	4000 x g
--------------------------	----------

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	28000
Height of case	125 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0298 cbm
Weight of case	4.03 kg
EAN inner box	4038917106962
EAN case	4038917011594

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT