



Product specification

Tube, 5 ml, (LxØ): 75 x 13 mm, PS



Product information

Product number	55.475.200
----------------	------------

Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 13 mm, PS
Product number	55.475.200
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PS, round base, transparent, push cap, suitable for Tecan Genesis, FE500, 490 piece(s)/StackPack

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	13 mm
Length of product	75 mm

Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	4000 x g

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / inner box	490
Piece(s) / case	1960
Piece(s) / pallet	109760
Height of inner box	193 mm
Width of inner box	77 mm
Depth of inner box	405 mm
Height of case	175 mm
Width of case	383 mm
Depth of case	428 mm
Volume of case	0.0287 cbm
Weight of case	6.8306 kg
EAN inner box	4038917158466
EAN case	4038917158459

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT