



Product specification

Tube, 5 ml, (LxØ): 75 x 12 mm, PS



Product information

Product number	55.476.013
----------------	------------

Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 12 mm, PS
Product number	55.476.013
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, natural, sterile, 25 piece(s)/bag

Material & colours

Color of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	natural

Size

Diameter	12 mm
Length of product	75 mm

Product characteristics

Base shape	Round base
Graduation	no
Cap	Ventilation cap
Closure type	Push cap
Centrifugation max (RCF)	4000 x g

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / inner box	25
Piece(s) / case	1000
Piece(s) / pallet	48000
Height of case	145 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0345 cbm
Weight of case	4.08 kg
EAN inner box	4038917107327
EAN case	4038917011785

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT