



Product specification

Tube, 5 ml, (LxØ): 75 x 12 mm, PS



Product information

Product number	55.476.061
----------------	------------

Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 12 mm, PS
Product number	55.476.061
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, suitable for Vitek and Vitek 2 from bioMérieux, 500 piece(s)/StackPack

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Size

Diameter	12 mm
Length of product	75 mm

Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	4000 x g

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	2000
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	120000
Height of inner box	176 mm
Width of inner box	77 mm
Depth of inner box	387 mm
Height of case	175 mm
Width of case	358 mm
Depth of case	398 mm
Volume of case	0.0249 cbm
Weight of case	6.2 kg
EAN inner box	4038917107631
EAN case	4038917011839

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT