



Product specification

Tube, 3.5 ml, (LxØ): 55 x 12 mm, PS



Product information

Product number	55.484.001
----------------	------------

Product descriptions

Name	Tube, 3.5 ml, (LxØ): 55 x 12 mm, PS
Product number	55.484.001
Product description	Tube, working volume: 3.5 ml, (LxØ): 55 x 12 mm, material: PS, round base, transparent, push cap, natural, cap assembled, sterile, 100 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent
Tube material	Polystyrene (PS)
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	natural

Size

Diameter	12 mm
Length of product	55 mm

Product characteristics

Base shape	Round base
Graduation	no
Cap	Standard cap

Closure type	Push cap
--------------	----------

Centrifugation max (RCF)	5000 x g
--------------------------	----------

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
---------------------------	--------------------------

Product category	no medical device IVD
------------------	-------------------------

Purity standard	sterile
-----------------	---------

Sterilisation	Electron irradiation
---------------	----------------------

Batched	yes
---------	-----

Packaging

Minimum order quantity	2000
------------------------	------

Piece(s) / inner box	100
----------------------	-----

Piece(s) / case	2000
-----------------	------

Piece(s) / pallet	64000
-------------------	-------

Height of case	145 mm
----------------	--------

Width of case	420 mm
---------------	--------

Depth of case	648 mm
---------------	--------

Volume of case	0.0395 cbm
----------------	------------

Weight of case	5.91 kg
----------------	---------

EAN inner box	4038917107846
---------------	---------------

EAN case	4038917012010
----------	---------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT