



## Product specification

**Tube, 5 ml, (LxØ): 75 x 13 mm, PP**



### Product information

Product number	55.525
----------------	--------

### Product descriptions

Name	Tube, 5 ml, (LxØ): 75 x 13 mm, PP
Product number	55.525
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, transparent, push cap, without cap, 500 piece(s)/StackPack

### Material & colours

Color of product	transparent
------------------	-------------

### Size

Diameter	13 mm
Length of product	75 mm

### Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	5000 x g

### Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	not regulated
CE certified	no CE mark
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	112000
Volume of case	0.0287 cbm
Weight of case	5.88 kg
EAN inner box	4038917082204
EAN case	4038917001298

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT