




Product specification

Tube, 11.5 ml, (LxØ): 100 x 15.7 mm, PP



Product information

Product number	55.539
----------------	--------

Product descriptions

Name	Tube, 11.5 ml, (LxØ): 100 x 15.7 mm, PP
Product number	55.539
Product description	Tube, working volume: 11.5 ml, (LxØ): 100 x 15.7 mm, material: PP, round base, transparent, push cap, 500 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	15.7 mm
Length of product	100 mm

Product characteristics

Base shape	Round base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	4000 x g

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Batched	no

Packaging

Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	42000
Height of inner box	60 mm
Width of inner box	535 mm
Depth of inner box	650 mm
Height of case	250 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0396 cbm
Weight of case	5.2 kg
EAN inner box	4038917082327
EAN case	4038917281997

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT