



Product specification

Tube, 12 ml, (LxØ): 110 x 17 mm, PS



Product information

Product number	57.462.001
----------------	------------

Product descriptions

Name	Tube, 12 ml, (LxØ): 110 x 17 mm, PS
Product number	57.462.001
Product description	Tube, working volume: 12 ml, (LxØ): 110 x 17 mm, material: PS, conical base, transparent, push cap, with graduations, sterile, 100 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent
Closure material	Low Density Polyethylene (LD-PE)

Size

Diameter	17 mm
Length of product	110 mm

Product characteristics

Base shape	Conical base
Graduation	yes
Closure type	Push cap
Centrifugation max (RCF)	2500 x g

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	24000
Height of case	145 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0345 cbm
Weight of case	4.35 kg
EAN inner box	4038917381048
EAN case	4038917381031

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT