




Product specification

Tube, 4.5 ml, (LxØ): 75 x 12 mm, PS



Product information

Product number	57.477
----------------	--------

Product descriptions

Name	Tube, 4.5 ml, (LxØ): 75 x 12 mm, PS
Product number	57.477
Product description	Tube, working volume: 4.5 ml, (LxØ): 75 x 12 mm, material: PS, conical base, transparent, push cap, without cap, 500 piece(s)/StackPack

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	12 mm
Length of product	75 mm

Product characteristics

Base shape	Conical base
Graduation	no
Closure type	Push cap
Centrifugation max (RCF)	4000 x g

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
---------------------------	--------------------------

Product category

In-vitro diagnostic (IVD)

CE certified

CE - manufacturer's self-declaration

Batched

yes

Packaging

Piece(s) / inner box

500

Piece(s) / case

2000

Piece(s) / pallet

120000

Height of inner box

176 mm

Width of inner box

77 mm

Depth of inner box

387 mm

Height of case

175 mm

Width of case

358 mm

Depth of case

398 mm

Volume of case

0.0249 cbm

Weight of case

5.84 kg

EAN inner box

4038917081856

EAN case

4038917000963

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT