



Product specification

Tube, 10 ml, (LxØ): 79 x 16 mm, PP



Product information

Product number	60.506
----------------	--------

Product descriptions

Name	Tube, 10 ml, (LxØ): 79 x 16 mm, PP
Product number	60.506
Product description	Tube, working volume: 10 ml, (LxØ): 79 x 16 mm, material: PP, round base, transparent, screw cap, natural, cap enclosed, 1,000 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	natural

Size

Diameter	16 mm
Length of product	79 mm

Product characteristics

Base shape	Round base
Graduation	no
Closure type	Screw cap
Centrifugation max (RCF)	4000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	1000
Piece(s) / case	1000
Piece(s) / pallet	36000
Height of case	359 mm
Width of case	400 mm
Depth of case	321 mm
Volume of case	0.046 cbm
Weight of case	4.96 kg
EAN inner box	4038917082037
EAN case	4038917035026

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT