



## Product specification

**Screw cap tube, 13 ml, (LxØ): 101 x 16.5 mm, PP**



### Product information

Product number	60.541.009
----------------	------------

### Product descriptions

Name	Screw cap tube, 13 ml, (LxØ): 101 x 16.5 mm, PP
Product number	60.541.009
Product description	Screw cap tube, working volume: 13 ml, (LxØ): 101 x 16.5 mm, material: PP, round base with skirted base, transparent, screw cap, without cap, 500 piece(s)/bag

### Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)

### Size

Diameter	16.5 mm
Length of product	101 mm

### Product characteristics

Base shape	Round base with skirted base
Graduation	no
Closure type	Screw cap
Autoclavable (max. °C)	121 °C

Centrifugation max (RCF)	4000 x g
--------------------------	----------

## Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	32000
Volume of case	0.0558 cbm
Weight of case	7.2 kg
EAN inner box	4038917110563
EAN case	4038917032438

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT