



## Product specification

**Screw cap tube, 10 ml, (LxØ): 79 x 16 mm, PP**



### Product information

Product number	60.551.001
----------------	------------

### Product descriptions

Name	Screw cap tube, 10 ml, (LxØ): 79 x 16 mm, PP
Product number	60.551.001
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 79 x 16 mm, material: PP, round base with skirted base, transparent, screw cap, natural, cap assembled, sterile, 100 piece(s)/bag

### Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	natural

### Size

Diameter	16 mm
Length of product	79 mm

### Product characteristics

Base shape	Round base with skirted base
Graduation	no

Cap	standard screw cap
Closure type	Screw cap
Centrifugation max (RCF)	4000 x g

## Purity & certification

Satisfies the requirement	IATA
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	36000
Height of case	464 mm
Width of case	257 mm
Depth of case	413 mm
Volume of case	0.0492 cbm
Weight of case	5.6 kg
EAN inner box	4038917112673
EAN case	4038917044103

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT