



Product specification

Screw cap tube, 120 ml, (LxØ): 114 x 44 mm, PP



Product information

Product number	60.597
----------------	--------

Product descriptions

Name	Screw cap tube, 120 ml, (LxØ): 114 x 44 mm, PP
Product number	60.597
Product description	Screw cap tube, working volume: 120 ml, (LxØ): 114 x 44 mm, material: PP, flat base, transparent, screw cap, natural, cap enclosed, 50 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	Polypropylene (PP)
Color of cap	natural

Size

Diameter	44 mm
Length of product	114 mm
Sample volume	120 ml

Product characteristics

Base shape	Flat base
------------	-----------

Graduation	no
Cap	standard screw cap
Closure type	Screw cap
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order quantity	250
Piece(s) / inner box	50
Piece(s) / case	250
Piece(s) / pallet	3000
Height of case	500 mm
Width of case	343 mm
Depth of case	563 mm
Volume of case	0.0966 cbm
Weight of case	6.48 kg
EAN inner box	4038917082785
EAN case	4038917001823

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT