



Product specification

Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm, conical false bottom, flat tube bottom, PP, without cap,



Product information

Product number	60.611.310
----------------	------------

Product descriptions

Name	Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm, conical false bottom, flat tube bottom, PP, without cap, 1,000 piece(s)/bag
Product number	60.611.310
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 92 x 15.3 mm, conical false bottom, flat tube bottom, transparent, material: PP, without cap, 1,000 piece(s)/bag, 1,000 piece(s)/carton

Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)

Size

Diameter	15.3 mm
Length excluding cap	92 mm
Length of product	92 mm
Height of product	92 mm

Product characteristics

Base shape	conical false bottom, flat tube bottom
Shape of inner tube	conical

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	1000
Piece(s) / pallet	30000
Volume of case	0.046 cbm
Weight of case	4.55 kg
EAN inner box	4038917559799
EAN case	4038917367905

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT