




## Product specification

**Screw cap tube, 2.5 ml, (LxØ): 75 x 13 mm, conical false bottom, rounded tube bottom, PP, without ca**



## Product information

Product number	60.614.065
----------------	------------

## Product descriptions

Name	Screw cap tube, 2.5 ml, (LxØ): 75 x 13 mm, conical false bottom, rounded tube bottom, PP, without cap, 100 piece(s)/bag
Product number	60.614.065
Product description	Screw cap tube, working volume: 2.5 ml, (LxØ): 75 x 13 mm, conical false bottom, rounded tube bottom, transparent, material: PP, with graduations, without cap, With injected scaling from 1.0 to 2.5 ml, 100 piece(s)/bag, 1,000 piece(s)/carton

## Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)

## Size

Diameter	13 mm
Length excluding cap	75 mm
Length of product	75 mm
Height of product	75 mm

## Product characteristics

Design	With injected scaling from 1.0 to 2.5 ml
Base shape	conical false bottom, rounded tube bottom
Graduation	yes

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	64000
Height of case	215 mm
Width of case	298 mm
Depth of case	398 mm
Volume of case	0.0255 cbm
Weight of case	2.82 kg
EAN inner box	4038917349765
EAN case	4038917349758

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT