



Product specification

Screw cap tube, 15 ml, (LxØ): 76 x 20 mm, PP



Product information

Product number	60.732
----------------	--------

Product descriptions

Name	Screw cap tube, 15 ml, (LxØ): 76 x 20 mm, PP
Product number	60.732
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 76 x 20 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap enclosed, 100 piece(s)/bag

Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	natural

Size

Diameter	20 mm
Length of product	76 mm

Product characteristics

Base shape	skirted conical base
Graduation	no
Cap	standard screw cap

Closure type	Screw cap
--------------	-----------

Centrifugation max (RCF)	4500 x g
--------------------------	----------

Purity & certification

Satisfies the requirement	ADR
---------------------------	-----

Product category	In-vitro diagnostic (IVD)
------------------	---------------------------

CE certified	CE - manufacturer's self-declaration
--------------	--------------------------------------

Batched	yes
---------	-----

Packaging

Piece(s) / inner box	100
----------------------	-----

Piece(s) / case	500
-----------------	-----

Piece(s) / pallet	18000
-------------------	-------

Height of case	265 mm
----------------	--------

Width of case	398 mm
---------------	--------

Depth of case	398 mm
---------------	--------

Volume of case	0.042 cbm
----------------	-----------

Weight of case	3.76 kg
----------------	---------

EAN inner box	4038917083805
---------------	---------------

EAN case	4038917002691
----------	---------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT