



## Product specification

**Screw cap tube, 15 ml, (LxØ): 76 x 20 mm, PP**



### Product information

Product number	60.732.001
----------------	------------

### Product descriptions

Name	Screw cap tube, 15 ml, (LxØ): 76 x 20 mm, PP
Product number	60.732.001
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 76 x 20 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap assembled, sterile, 100 piece(s)/bag

### Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	natural

### Size

Diameter	20 mm
Length of product	76 mm

### Product characteristics

Base shape	skirted conical base
Graduation	no
Cap	standard screw cap

Closure type	Screw cap
Centrifugation max (RCF)	4500 x g

## Purity & certification

Satisfies the requirement	ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	20000
Height of case	165 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0393 cbm
Weight of case	3.78 kg
EAN inner box	4038917119603
EAN case	4038917016698

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT