



Product specification

Screw cap tube, 10 ml, (LxØ): 97 x 16 mm, PS



Product information

Product number	60.9921.829
----------------	-------------

Product descriptions

Name	Screw cap tube, 10 ml, (LxØ): 97 x 16 mm, PS
Product number	60.9921.829
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 97 x 16 mm, material: PS, skirted conical base, transparent, screw cap, yellow, cap assembled, sterile, 100 piece(s)/bag

Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	yellow

Size

Diameter	16 mm
Length of product	97 mm

Product characteristics

Base shape	skirted conical base
Graduation	no
Closure type	Screw cap

Centrifugation max (RCF)

4000 x g

Purity & certification

Satisfies the requirement	ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	30000
Height of inner box	50 mm
Width of inner box	264 mm
Depth of inner box	410 mm
Height of case	264 mm
Width of case	324 mm
Depth of case	608 mm
Volume of case	0.052 cbm
Weight of case	8.5 kg
EAN inner box	4038917151122
EAN case	4038917027700

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT