



Product specification

Screw cap tube, 12 ml, (LxØ): 99 x 16 mm, PC



Product information

| | |
|----------------|-------------|
| Product number | 60.9922.936 |
|----------------|-------------|

Product descriptions

| | |
|---------------------|---|
| Name | Screw cap tube, 12 ml, (LxØ): 99 x 16 mm, PC |
| Product number | 60.9922.936 |
| Product description | Screw cap tube, working volume: 12 ml, (LxØ): 99 x 16 mm, material: PC, round base, transparent, screw cap, white, cap enclosed, 1,000 piece(s)/bag |

Material & colours

| | |
|------------------|--------------------|
| Color of product | transparent |
| Tube material | Polycarbonate (PC) |
| Closure material | Polypropylene (PP) |
| Color of cap | white |

Size

| | |
|-------------------|-------|
| Diameter | 16 mm |
| Length of product | 99 mm |

Product characteristics

| | |
|------------|------------|
| Base shape | Round base |
| Graduation | no |

| | |
|--------------------------|--------------------|
| Cap | standard screw cap |
| Closure type | Screw cap |
| Autoclavable (max. °C) | 121 °C |
| Centrifugation max (RCF) | 8800 x g |

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | ADR |
| Product category | no medical device IVD |
| Latex-free | yes |
| Batched | yes |
| Safety data sheet | no |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 1000 |
| Piece(s) / case | 1000 |
| Piece(s) / pallet | 24000 |
| Height of case | 264 mm |
| Width of case | 373 mm |
| Depth of case | 604 mm |
| Volume of case | 0.0595 cbm |
| Weight of case | 8.4 kg |
| EAN inner box | 4038917153423 |
| EAN case | 4038917027816 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT