



## Product specification

**Screw cap tube, 12 ml, (LxØ): 99 x 16 mm, PC**



### Product information

Product number	60.9922.937
----------------	-------------

### Product descriptions

Name	Screw cap tube, 12 ml, (LxØ): 99 x 16 mm, PC
Product number	60.9922.937
Product description	Screw cap tube, working volume: 12 ml, (LxØ): 99 x 16 mm, material: PC, round base, transparent, screw cap, yellow, cap assembled, sterile, 100 piece(s)/bag

### Material & colours

Product material	Polycarbonate (PC)
Color of product	transparent
Closure material	Polypropylene (PP)
Color of cap	yellow

### Size

Diameter	16 mm
Length of product	99 mm

### Product characteristics

Base shape	Round base
Graduation	no

Closure type	Screw cap
Centrifugation max (RCF)	8800 x g

## Purity & certification

Satisfies the requirement	ADR
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	30000
Height of case	264 mm
Width of case	373 mm
Depth of case	604 mm
Volume of case	0.0595 cbm
Weight of case	8.4 kg
EAN inner box	4038917153430
EAN case	4038917027823

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT