



## Product specification

**Screw cap tube, 8 ml, (LxØ): 94 x 14 mm, PC**



### Product information

Product number	60.9922.941
----------------	-------------

### Product descriptions

Name	Screw cap tube, 8 ml, (LxØ): 94 x 14 mm, PC
Product number	60.9922.941
Product description	Screw cap tube, working volume: 8 ml, (LxØ): 94 x 14 mm, material: PC, round base, transparent, screw cap, yellow, cap assembled, sterile, 100 piece(s)/bag

### Material & colours

Color of product	transparent
Closure material	Polypropylene (PP)
Color of cap	yellow

### Size

Diameter	14 mm
Length of product	94 mm

### Product characteristics

Base shape	Round base
Graduation	no
Closure type	Screw cap

Centrifugation max (RCF)

8800 x g

## Purity & certification

Satisfies the requirement	ADR
Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	36000
Height of case	264 mm
Width of case	324 mm
Depth of case	608 mm
Volume of case	0.052 cbm
Weight of case	7.2 kg
EAN inner box	4038917153478
EAN case	4038917027847

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT