



## Product specification

**Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP,  
cap enclosed, 1,000 piece(s)/bag**



## Product information

Product number	62.504.038
----------------	------------

## Product descriptions

Name	Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, cap enclosed, 1,000 piece(s)/bag
Product number	62.504.038
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, brown, material: PP, with print, label/print: black, with graduations, cap enclosed, brown, 1,000 piece(s)/bag, 1,000 piece(s)/carton

## Material & colours

Color of product	brown
Tube material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	brown

## Size

Diameter	13 mm
Length of product	75 mm

## Product characteristics

Base shape	Round base
Label/Print	With print
Color of print/label	black
Graduation	yes
Light protection	yes
Closure type	Screw cap

## Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Piece(s) / inner box	1000
Piece(s) / case	1000
Piece(s) / pallet	56000
Height of case	120 mm
Width of case	433 mm
Depth of case	568 mm
Volume of case	0.0295 cbm
Weight of case	4.39 kg
EAN inner box	4038917412971
EAN case	4038917412964

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT