



Product specification

**Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP,
cap enclosed, 100 piece(s)/bag**



Product information

Product number	62.504.039
----------------	------------

Product descriptions

Name	Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, cap enclosed, 100 piece(s)/bag
Product number	62.504.039
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, brown, material: PP, with print, label/print: black, with graduations, cap enclosed, brown, 100 piece(s)/bag, 1,000 piece(s)/carton

Material & colours

Color of product	brown
Tube material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	brown

Size

Diameter	13 mm
Length of product	75 mm

Product characteristics

Base shape	Round base
Label/Print	With print
Color of print/label	black
Graduation	yes
Light protection	yes
Closure type	Screw cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	40000
Height of case	205 mm
Width of case	408 mm
Depth of case	458 mm
Volume of case	0.0383 cbm
Weight of case	4.39 kg
EAN inner box	4038917413008
EAN case	4038917412995

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT