



## Product specification

**Tube, 13 ml, (LxØ): 100 x 16 mm, PP, with print**



### Product information

Product number	62.515.006
----------------	------------

### Product descriptions

Name	Tube, 13 ml, (LxØ): 100 x 16 mm, PP, with print
Product number	62.515.006
Product description	Tube, working volume: 13 ml, (LxØ): 100 x 16 mm, material: PP, round base, transparent, push cap, natural, cap assembled, with print, label/print: white, with graduations, sterile, 25 piece(s)/bag

### Material & colours

Color of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	natural

### Size

Diameter	16 mm
Length of product	100 mm

### Product characteristics

Base shape	Round base
Label/Print	With print
Color of print/label	white

Graduation	yes
Cap	Ventilation cap
Closure type	Push cap
Centrifugation max (RCF)	3000 x g

## Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Piece(s) / inner box	25
Piece(s) / case	500
Piece(s) / pallet	24000
Height of case	145 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0345 cbm
Weight of case	3.05 kg
EAN inner box	4038917108836
EAN case	4038917012263

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT