



## Product specification

**Screw cap tube, 30 ml, (LxØ): 107 x 25 mm, PP, with print**



## Product information

Product number	62.543.001
----------------	------------

## Product descriptions

Name	Screw cap tube, 30 ml, (LxØ): 107 x 25 mm, PP, with print
Product number	62.543.001
Product description	Screw cap tube, working volume: 30 ml, (LxØ): 107 x 25 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap assembled, with print, label/print: black, with graduations, sterile, 50 piece(s)/bag

## Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	natural

## Size

Diameter	25 mm
Length of product	107 mm

## Product characteristics

Base shape	skirted conical base
Label/Print	With print
Color of print/label	black

Graduation	yes
Cap	standard screw cap
Closure type	Screw cap
Centrifugation max (RCF)	8800 x g

## Purity & certification

Satisfies the requirement	ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Piece(s) / inner box	50
Piece(s) / case	400
Piece(s) / pallet	11200
Height of case	255 mm
Width of case	400 mm
Depth of case	590 mm
Volume of case	0.0602 cbm
Weight of case	5.59 kg
EAN inner box	4038917111638
EAN case	4038917013284

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT