



Product specification

Screw cap tube, 3.5 ml, (LxØ): 66 x 11.5 mm, PP, with print



Product information

Product number	62.549.004
----------------	------------

Product descriptions

Name	Screw cap tube, 3.5 ml, (LxØ): 66 x 11.5 mm, PP, with print
Product number	62.549.004
Product description	Screw cap tube, working volume: 3.5 ml, (LxØ): 66 x 11.5 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap enclosed, with print, label/print: white, with graduations, 1,000 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	natural

Size

Diameter	11.5 mm
Length of product	66 mm

Product characteristics

Base shape	skirted conical base
Label/Print	With print

Color of print/label	white
Graduation	yes
Closure type	Screw cap
Centrifugation max (RCF)	4000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Batched	no

Packaging

Minimum order quantity	50000
Piece(s) / inner box	1000
Piece(s) / case	2000
Piece(s) / pallet	84000
Height of case	250 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0396 cbm
Weight of case	5.7 kg
EAN inner box	4038917112024
EAN case	4038917013468

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT