



## Product specification

**Screw cap tube, 30 ml, (LxØ): 84 x 30 mm, PP, with print**



## Product information

Product number	62.555
----------------	--------

## Product descriptions

Name	Screw cap tube, 30 ml, (LxØ): 84 x 30 mm, PP, with print
Product number	62.555
Product description	Screw cap tube, working volume: 30 ml, (LxØ): 84 x 30 mm, material: PP, flat base, transparent, screw cap, natural, cap enclosed, with print, label/print: white, with graduations, 50 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	Polypropylene (PP)
Color of cap	natural

## Size

Diameter	30 mm
Length of product	84 mm

## Product characteristics

Base shape	Flat base
Label/Print	With print

Color of print/label	white
Graduation	yes
Closure type	Screw cap
Autoclavable (max. °C)	121 °C

## Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

## Packaging

Minimum order quantity	450
Piece(s) / inner box	50
Piece(s) / case	450
Piece(s) / pallet	10800
Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	5.8302 kg
EAN inner box	4038917082433
EAN case	4038917079754

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT