




Product specification

Screw cap tube, 10 ml, (LxØ): 92 x 15.3 mm, PP, with print



Product information

Product number	62.610.018
----------------	------------

Product descriptions

Name	Screw cap tube, 10 ml, (LxØ): 92 x 15.3 mm, PP, with print
Product number	62.610.018
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 92 x 15.3 mm, material: PP, round base, transparent, screw cap, natural, cap assembled, with print, label/print: green, "CSF/Liquor", with graduations, sterile, 1 piece(s)/blister

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	natural

Size

Diameter	15.3 mm
Length of product	92 mm

Product characteristics

Base shape	Round base
------------	------------

Label/Print	With print
Color of print/label	green
Graduation	yes
Closure type	Screw cap
Centrifugation max (RCF)	3000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	13500
Height of inner box	171 mm
Width of inner box	151 mm
Depth of inner box	351 mm
Height of case	195 mm
Width of case	368 mm
Depth of case	778 mm
Volume of case	0.0558 cbm
Weight of case	4.23 kg
EAN inner box	4038917233132
EAN case	4038917233101

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



SARSTEDT

This document was prepared by EDP support and is valid without signature.