



Product specification

Screw cap tube, 3.5 ml, (LxØ): 92 x 13 mm, conical false bottom, flat tube bottom, PP, cap assembled



Product information

| | |
|----------------|--------|
| Product number | 62.613 |
|----------------|--------|

Product descriptions

| | |
|---------------------|--|
| Name | Screw cap tube, 3.5 ml, (LxØ): 92 x 13 mm, conical false bottom, flat tube bottom, PP, cap assembled, 100 piece(s)/bag |
| Product number | 62.613 |
| Product description | Screw cap tube, working volume: 3.5 ml, (LxØ): 92 x 13 mm, conical false bottom, flat tube bottom, transparent, material: PP, with graduations, cap assembled, 100 piece(s)/bag, 1,000 piece(s)/carton |

Material & colours

| | |
|------------------|-------------|
| Color of product | transparent |
|------------------|-------------|

Size

| | |
|-------------------|-------|
| Diameter | 13 mm |
| Length of product | 92 mm |

Product characteristics

| | |
|------------|--|
| Base shape | conical false bottom, flat tube bottom |
|------------|--|

Graduation

yes

Purity & certification

| | |
|---------------------------|---------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 100 |
| Piece(s) / case | 1000 |
| Piece(s) / pallet | 20000 |
| Height of inner box | 20 mm |
| Width of inner box | 330 mm |
| Depth of inner box | 340 mm |
| Height of case | 210 mm |
| Width of case | 410 mm |
| Depth of case | 460 mm |
| Volume of case | 0.0396 cbm |
| Weight of case | 4.34 kg |
| EAN inner box | 4038917082891 |
| EAN case | 4038917079136 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT