




Product specification

Tube, 10 ml, (LxØ): 100 x 16 mm, PP, with print



Product information

Product number	62.9924.272
----------------	-------------

Product descriptions

Name	Tube, 10 ml, (LxØ): 100 x 16 mm, PP, with print
Product number	62.9924.272
Product description	Tube, working volume: 10 ml, (LxØ): 100 x 16 mm, material: PP, conical base, transparent, push cap, with print, label/print: white, with graduations, 100 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	16 mm
Length of product	100 mm

Product characteristics

Base shape	Conical base
Label/Print	With print
Color of print/label	white
Graduation	yes
Closure type	Push cap

Centrifugation max (RCF)

4500 x g

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	30000
Height of case	320 mm
Width of case	260 mm
Depth of case	608 mm
Volume of case	0.0506 cbm
Weight of case	5 kg
EAN inner box	4038917155120
EAN case	4038917037525

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT