




## Product specification

**Tube, 10 ml, (LxØ): 100 x 16 mm, PP, with print**



## Product information

Product number	62.9924.272
----------------	-------------

## Product descriptions

Name	Tube, 10 ml, (LxØ): 100 x 16 mm, PP, with print
------	---

Product number	62.9924.272
----------------	-------------

Product description	Tube, working volume: 10 ml, (LxØ): 100 x 16 mm, material: PP, conical base, transparent, push cap, with print, label/print: white, with graduations, 100 piece(s)/bag
---------------------	--

## Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

Color of product	transparent
------------------	-------------

## Size

Diameter	16 mm
----------	-------

Length of product	100 mm
-------------------	--------

## Product characteristics

Base shape	Conical base
------------	--------------

Label/Print	With print
-------------	------------

Color of print/label	white
----------------------	-------

Graduation	yes
------------	-----

Closure type	Push cap
Centrifugation max (RCF)	4500 x g

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	30000
Height of case	264 mm
Width of case	324 mm
Depth of case	608 mm
Volume of case	0.052 cbm
Weight of case	5 kg
EAN inner box	4038917155120
EAN case	4038917037525

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT