




## Product specification

**Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP, with paper label**



## Product information

Product number	63.9922.252
----------------	-------------

## Product descriptions

Name	Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP, with paper label
Product number	63.9922.252
Product description	Screw cap tube, working volume: 25 ml, (LxØ): 90 x 25 mm, material: PP, skirted conical base, transparent, screw cap, white, cap assembled, with paper label, label/print: white/black, Patient information, with graduations, 500 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	white

## Size

Diameter	25 mm
Length of product	90 mm

## Product characteristics

Base shape	skirted conical base
Label/Print	With paper label
Color of print/label	white/black
Graduation	yes
Closure type	Screw cap
Centrifugation max (RCF)	4500 x g

## Purity & certification

Satisfies the requirement	IATA
Product category	diagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order quantity	500
Piece(s) / inner box	500
Piece(s) / case	500
Piece(s) / pallet	15000
Height of case	264 mm
Width of case	324 mm
Depth of case	608 mm
Volume of case	0.052 cbm
Weight of case	4.98 kg
EAN inner box	4038917151559
EAN case	4038917027779

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT