



# Product specification

Cap, yellow, suitable for tubes Ø 12-17 mm



## Product information

Product number	65.1379.004
----------------	-------------

## Product descriptions

Name	Cap, yellow, suitable for tubes Ø 12-17 mm
Product number	65.1379.004
Product description	Cap, yellow, suitable for tubes Ø 12-17 mm, 1,000 piece(s)/bag

## Material & colours

Product material	Low Density Polyethylene (LD-PE)
Color of product	yellow

## Product characteristics

Cap	Standard cap
-----	--------------

## Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Minimum order quantity	10000
Piece(s) / inner box	1000
Piece(s) / case	10000

Piece(s) / pallet	240000
Height of case	464 mm
Width of case	257 mm
Depth of case	413 mm
Volume of case	0.0492 cbm
Weight of case	5.36 kg
EAN inner box	4038917127998
EAN case	4038917021524

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT