



Product specification

Cap, green, suitable for tubes Ø 10-17 mm



Product information

Product number	65.1379.005
----------------	-------------

Product descriptions

Name	Cap, green, suitable for tubes Ø 10-17 mm
Product number	65.1379.005
Product description	Cap, green, suitable for tubes Ø 10-17 mm, 1,000 piece(s)/bag

Material & colours

Color of product	green
------------------	-------

Product characteristics

Cap	Standard cap
Closure type	Push cap

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	360000

Volume of case	0.0492 cbm
Weight of case	5.35 kg
EAN inner box	4038917128001
EAN case	4038917021531

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT