




Product specification

Screw cap, red, suitable for tubes Ø 16-16.5 mm



Product information

Product number	65.176.001
----------------	------------

Product descriptions

Name	Screw cap, red, suitable for Tubes Ø 16-16.5 mm
Product number	65.176.001
Product description	Screw cap, red, suitable for tubes Ø 16-16.5 mm, 1,001 piece(s)/bag

Material & colours

Color of product	red
------------------	-----

Product characteristics

Closure type	Screw cap
--------------	-----------

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	5000
Piece(s) / pallet	200000
Volume of case	0.0415 cbm

Weight of case	7.15 kg
EAN inner box	4038917318044
EAN case	4038917318037

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT