



Product specification

Screw cap, purple, suitable for tubes 82 x 13 mm



Product information

Product number	65.177.042
----------------	------------

Product descriptions

Name	Screw cap, purple, suitable for Tubes 82 x 13 mm
Product number	65.177.042
Product description	Screw cap, purple, suitable for tubes 82 x 13 mm, 1,000 piece(s)/bag

Material & colours

Color of product	purple
------------------	--------

Product characteristics

Closure type	Screw cap
--------------	-----------

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	480000
Volume of case	0.0366 cbm

Weight of case	7.76 kg
EAN inner box	4038917101707
EAN case	4038917028240

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT