




Product specification

Colour-coded inserts, for CryoPure tubes, white



Product information

Product number	65.386
----------------	--------

Product descriptions

Name	Colour-coded inserts, for CryoPure tubes, white
Product number	65.386
Product description	Colour-coded inserts, for CryoPure tubes, white, 100 piece(s)/bag

Material & colours

Color of product	white
------------------	-------

Product characteristics

Graduation	no
------------	----

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	3000
Piece(s) / pallet	972000

Height of inner box	10 mm
Width of inner box	125 mm
Depth of inner box	175 mm
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	0.681 kg
EAN inner box	4038917081429
EAN case	4038917077767

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT