


# Product specification

## Colour-coded inserts, for CryoPure tubes, green



### Product information

Order number	65.386.005
--------------	------------

### Product descriptions

Name	Colour-coded inserts, for CryoPure tubes, green
Product number	65.386.005
Product description	Colour-coded inserts, for CryoPure tubes, green, 100 piece(s)/bag

### Product characteristics

Graduation	no
Suitable for	CryoPure tubes

### Material & colours

Product material	Polystyrene (PS)
Color of product	green

### Purity & certification

Accordance with MDD or IVDD	no medical device   IVD
Lot or serial number	yes

### Packaging

Piece(s) / inner box	100
Piece(s) / case	3000
Piece(s) / pallet	972000
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	0.75 kg

**This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.**

**This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.**

This document was prepared by EDP support and is valid without signature.

