




## Product specification

### Colour-coded inserts, for CryoPure tubes, blue



## Product information

Product number	65.386.006
----------------	------------

## Product descriptions

Name	Colour-coded inserts, for CryoPure tubes, blue
Product number	65.386.006
Product description	Colour-coded inserts, for CryoPure tubes, blue, 100 piece(s)/bag

## Material & colours

Product material	Polystyrene (PS)
Color of product	blue

## Product characteristics

Graduation	no
------------	----

## Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Minimum order quantity	100
Piece(s) / inner box	100
Piece(s) / case	3000

Piece(s) / pallet	972000
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	0.681 kg
EAN inner box	4038917274623
EAN case	4038917274616

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT