

Product specification

Cover for Pressure infusion cuff FR



Product information

Order number	97.8100.210
--------------	-------------

Product descriptions

Name	Cover for Pressure infusion cuff FR
------	-------------------------------------

Product number	97.8100.210
----------------	-------------

Product description	Cover, Replacement part for FR pressure infusion cuff
---------------------	---

Product characteristics

Suitable for	Pressure infusion cuff
--------------	------------------------

Size

Width of product	14.5 cm
------------------	---------

Length of product	27 cm
-------------------	-------

Height of product	2 cm
-------------------	------

Max. volume	1000 ml
-------------	---------

Material & colours

Color of product	red
------------------	-----

Purity & certification

Product category	no medical device IVD
------------------	-------------------------

Batched	no
---------	----

Packaging

Volume of case	0.0039 cbm
----------------	------------

Weight of case	0.97 kg
----------------	---------

EAN inner box	4028339311228
---------------	---------------

EAN case	4028339311211
----------	---------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

