




## Product specification

**Push cap, orange, suitable for tubes Ø 16-17 mm**



### Product information

Product number	65.793.516
----------------	------------

### Product descriptions

Name	Push cap, orange, suitable for Tubes Ø 16-17 mm
Product number	65.793.516
Product description	Push cap, orange, suitable for tubes Ø 16-17 mm, 1,000 piece(s)/bag

### Material & colours

Color of product	orange
------------------	--------

### Product characteristics

Closure type	Push cap
--------------	----------

### Purity & certification

Product category	no medical device   IVD
Batched	yes

### Packaging

Piece(s) / inner box	1000
Piece(s) / case	5000
Piece(s) / pallet	200000
Volume of case	0.0414 cbm

Weight of case	7.62 kg
EAN inner box	4038917120869
EAN case	4038917033718

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT