

# Product specification

## Push cap, yellow, suitable for tubes Ø 13 mm



### Product information

Product number	65.806.307
----------------	------------

### Product descriptions

Name	Push cap, yellow, suitable for Tubes Ø 13 mm
Product number	65.806.307
Product description	Push cap, yellow, suitable for tubes Ø 13 mm, 1,000 piece(s)/bag

### Product characteristics

Type of cap	Push cap
Suitable for	Tubes Ø 13 mm

### Material & colours

Product material	Low Density Polyethylene (LD-PE)
Color of product	yellow

### Purity & certification

Product category	no medical device   IVD
Batched	yes

### Packaging

Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	400000
Volume of case	0.0382 cbm
Weight of case	6.98 kg
EAN inner box	4038917368872
EAN case	4038917368865

**This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.**

**This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.**

This document was prepared by EDP support and is valid without signature.

