



Product specification

Push cap, yellow, suitable for tubes Ø 12 mm



Product information

Product number	65.809.301
----------------	------------

Product descriptions

Name	Push cap, yellow, suitable for Tubes Ø 12 mm
Product number	65.809.301
Product description	Push cap, yellow, suitable for tubes Ø 12 mm, 1,000 piece(s)/bag

Material & colours

Color of product	yellow
------------------	--------

Product characteristics

Closure type	Push cap
--------------	----------

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	400000
Volume of case	0.0382 cbm

Weight of case	6.81 kg
EAN inner box	4038917369145
EAN case	4038917369138

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT