

Product specification

Push cap, red, suitable for tubes Ø 12 mm



Product information

Order number	65.809.306
--------------	------------

Product descriptions

Name	Push cap, red, suitable for Tubes Ø 12 mm
Product number	65.809.306
Product description	Push cap, red, suitable for tubes Ø 12 mm, 1000 piece(s)/bag

Product characteristics

Type of cap	Push cap
Suitable for	Tubes Ø 12 mm

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Color of product	red

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	400000
Volume of case	0.0382 cbm
Weight of case	6.82 kg
EAN inner box	4038917369206
EAN case	4038917369190

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

