



## Product specification

**Quick-release cap 1.60–2.10 mm for blood gas capillaries**

**100 µl / 125 µl, white**



## Product information

Product number	65.935.205
----------------	------------

## Product descriptions

Name	Quick-release cap 1.60–2.10 mm for blood gas capillaries 100 µl / 125 µl, white
Product number	65.935.205
Product description	Quick-release cap 1.60–2.10 mm for blood gas capillaries 100 µl / 125 µl, white, blood gas analysis, 500 piece(s)/bag

## Material & colours

Color of product	white
------------------	-------

## Product characteristics

Type of collection	capillary
--------------------	-----------

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	32000
Piece(s) / pallet	2560000
Height of case	180 mm
Width of case	288 mm
Depth of case	378 mm
Volume of case	0.0196 cbm
Weight of case	8.416 kg
EAN inner box	4038917234559
EAN case	4038917234542

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT