




Product specification

PCR strip of lids, Biosphere® plus, PP, highly transparent, flat cap



Product information

Product number	65.989
----------------	--------

Product descriptions

Name	PCR strip of lids, Biosphere® plus, PP, highly transparent, flat cap
Product number	65.989
Product description	PCR strip of lids, for PCR plates (REF 72.1978/.010, 72.1979/.010/.003/.201/.700, 72.1980/.010/.201/.600/.700, 72.1981/.010) & PCR strips (REF 72.985.002/.092/.992), Biosphere® plus, material: PP, highly transparent, flat cap, suitable for qPCR, 12 piece(s)/bag

Material & colours

Color of product	highly transparent
------------------	--------------------

Size

Diameter	4.47 mm
Width of product	7.3 mm
Length of product	77.2 mm
Height of product	3.25 mm

Product characteristics

Cap	flat cap
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Pierceable	no
Suitable for qPCR	yes

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	12
Piece(s) / case	240
Piece(s) / pallet	77760
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	0.2599 kg
EAN inner box	4038917281485
EAN case	4038917281478

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT