


Product specification

Filter tip, 2.5 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/box



Product information

Product number	70.1130.217
----------------	-------------

Product descriptions

Name	Filter tip, 2.5 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/Box
Product number	70.1130.217
Product description	Filter tip, working volume: 2.5 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for Eppendorf, Gilson and Finnpiquette and products of identical design, 96 piece(s)/box

Product characteristics

Filter	yes
Graduation	yes
Suitable for	Eppendorf, Gilson and Finnpiquette and products of identical design

Size

Length of product	31 mm
Max. volume	2.5 µl

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Filter material	Polyethylene (PE)

Purity & certification

Product category	no medical device IVD
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	157440
Height of inner box	131 mm
Width of inner box	89 mm
Depth of inner box	314 mm
Height of case	150 mm
Width of case	326 mm
Depth of case	370 mm
Volume of case	0.0181 cbm
Weight of case	2.8493 kg
EAN inner box	4038917327282
EAN case	4038917327275

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

