




## Product specification

**Pipette tip, 5 ml, transparent, 250 piece(s)/bag**



### Product information

Product number	70.1183.001
----------------	-------------

### Product descriptions

Name	Pipette tip, 5 ml, transparent, 250 piece(s)/bag
Product number	70.1183.001
Product description	Pipette tip, working volume: 5 ml, transparent, filling level rings, suitable for Gilson, 250 piece(s)/bag

### Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

### Size

Length of product	123 mm
Max. volume	5 ml

### Product characteristics

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

### Purity & certification

Product category	no medical device   IVD
------------------	-------------------------

Batched

yes

## Packaging

Minimum order quantity	2500
Piece(s) / inner box	250
Piece(s) / case	2500
Piece(s) / pallet	25000
Height of case	698 mm
Width of case	464 mm
Depth of case	298 mm
Volume of case	0.0966 cbm
Weight of case	7.3 kg
EAN inner box	4038917126441
EAN case	4038917574587

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT