




Product specification

Pipette tip, 5 ml, transparent, 250 piece(s)/bag



Product information

Product number	70.1183.002
----------------	-------------

Product descriptions

Name	Pipette tip, 5 ml, transparent, 250 piece(s)/bag
Product number	70.1183.002
Product description	Pipette tip, working volume: 5 ml, transparent, filling level rings, suitable for Eppendorf, 250 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Length of product	123 mm
Max. volume	5 ml

Product characteristics

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	not regulated
CE certified	no CE mark

Batched	yes
---------	-----

Packaging

Piece(s) / inner box	250
Piece(s) / case	2500
Piece(s) / pallet	25000
Height of inner box	65 mm
Width of inner box	325 mm
Depth of inner box	520 mm
Height of case	450 mm
Width of case	300 mm
Depth of case	500 mm
Volume of case	0.0675 cbm
Weight of case	7.3 kg
EAN inner box	4038917126458
EAN case	4038917574594

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT